



Prior to analyses of Phosphorylated Proteins & neurotransmitters !

- All control of the instrument is through the Touch Panel Display on the front panel .
- Power level setting has become easier and more reproducible due to a digital setting
- Switching power supply has been employed to comply with line voltage fluctuation
- When ordering, please specify either 3-Phase 190-240V 30 amp or 380-440V 20 amp
- As air-cooling method is employed, water is no longer required for cooling.
- Muromachi Microwave Fixation Systems are safely designed so that the microwave leakage will not exceed $1\text{mW}/\text{cm}^2$