



LC 80, as floor loading cart and sterile containers washer disinfectant, has been developed for any Health Care application where large items are required to be cleaned and thermo disinfected. It can treat cleaning theatre transport and storage trolleys, sterile containers, operating tables, bed frames etc.

- The flat chamber floor provides a safe surface to walk during loading-unloading phase. Multiple oscillating pipes grants excellent washing results and the incorporated tilting rails assure optimum washing-drying and items unloading.


- LC 80 fulfill the European EN ISO15883-1/2 and UK HTM 2030 requirements.



-Flush self cleaning filter for the re-circulation of the wash water. Optional.



-Multiple oscillating pipes with numerous nozzles grant excellent washing results.




-Multiple rails for washing cycle and unloading of trolley.



- Door made in double HST tempered glass.
- Stainless steel door with glass window available as option.



-The STEELCOTRONIC system enables you to optimize the washing process and to personalize the thermo disinfection programmes built into the machine. In addition it allows the user to trace the principal validated parameters according to the European rules regarding traceability of contaminated items. The control screen allows viewing of all programme phases throughout the machines function and also highlights residual time.



-Patented manifold system for trolley docking with spray arms available - Optional.