



MK-400 has been developed for the auditory response test of the small animals. An animal is placed in a cage and given a sound established by the Audio Generator. The parameters are observed by the operator.

### Features

- Depress the START button after setting frequency, volume and sound duration, and the sound pressure will be produced accordingly.
- Sound can be set in three (3) levels ----- High (H), Middle (M), Low (L) and frequency of each level is determined at the time of manufacturing under the directions of the researcher.
- Sound duration can be set for the range of 0.1 second - 9.9 seconds by 0.1 sec steps

### Specification

Model	MK-400
AG-400 Audio Generator	1 (200W x 230D x 90mmH)
Speaker	1 (220W x 125D x 220mmH)
RC-400 Rat Cage	1 (230W x 205D x 230mmH)