

LAB 500C

FEEDING BOTTLES WASHERS



Washing injection System and air drying on two levels



Chamber volume **171 Lt**
Basket volume **151 Lt**

- LAB 500 C is specially designed for laboratories with limited space. A powerful filtering air drying system is a feature of this under bench laboratory glassware washer. Washing and air drying on 2 levels, normally seen on large machines, has been incorporated into this small model.

LAB 600

FEEDING BOTTLES WASHERS



Washing injection System and air drying on Three levels



Chamber volume **200 Lt**
Basket volume **170 Lt**


- LAB 600 front loading laboratory glassware features are the adjustable upper washing/drying levels and a powerful built-in HEPA filtered forced air drying system.
- The upper level can be placed at two different water/air connections depending on the height of the loaded glassware.

LAB 610

FEEDING BOTTLES WASHERS



Washing injection System and air drying on Four levels



Chamber volume **250 Lt**
Basket volume **220 Lt**

- LAB 610 glassware is a larger model, frontal loading, with HEPA-filtered forced air drying, four washing/drying positions and three washing/drying levels.
- The upper level can be placed at two different water/air connections depending on the height of the loaded glassware. This machine can be used also to clean large bottles (20,50,60 Lt.).