



Features

- Step-Through Test System has been designed to test animal memory.
- The test cage is divided into the two compartments, the light and the dark by the guillotine door.
- The animal is placed in the light compartment and allowed to explore there for a preset period of time so that it may be familiar with its environment.
- Then the door is opened and it tries to get in the dark room since it prefers a dark place. As soon as it gets into the dark, an electrical shock is given manually.

Specification

Model		MST-01SM FOR MOUSE	MST-01SR FOR RATS
Dimensions	Dark Room	140W x 175D x 150 mmH	280W x 295D x 300 mmH
	Light Room	90W x 115D x 150 mmH	240W x 245D x 300 mmH
	Overall	360W x 180D x 175 mmH	610W x 320D x 340 mmH
Grid	Diameter	4mm O.D.	5 mm O.D.
	Pitch To Pitch	10 mm	15 mm
	Grid To Grid	6 mm	10 mm
Weight		Approx. 2 kg	Approx. 8 kg